




5000W STRONG POWER
4 HANDLES WORK SIMULTANEOUSLY

+25% Muscle MORE

-30% Fat LESS


COMBINED 30-MIN. TREATMENT



RF

High Intensity Focused Electro-Magnetic

-30% Fat LESS



+25% Muscle MORE





EMS machine is based on an applicator simultaneously emitting synchronized RF and High Intensity Focused Electro-Magnetic energies.


Due to the radiofrequency heating, the muscle temperature quickly raises by several degrees. This prepares muscles for exposure to stress, similar to what a warm up activity does before any workout. In less than 4 minutes, the temperature in subcutaneous fat reaches levels that cause apoptosis, i.e. fat cells are permanently damaged and slowly removed from the body. Clinical studies showed on average a 30% reduction in subcutaneous fat.*

Bypassing the brain limitations, High Intensity Focused Electro-Magnetic energy contracts the muscle fibers in the area at intensities that are not achievable during voluntary workout. Extreme stress forces the muscle to adapt, resulting in an increase in the number and growth of muscle fibers and cells. Clinical studies showed on average a 25% muscle volume increase.*


Mechanism of muscle contraction



Nerve impulses set off a biochemical reaction that causes myosin to stick to actin



EMS machine induces very strong (supramaximal) muscle contractions, not achievable through workout. When exposed to these strong contractions, the muscle tissue is forced to adapt to such extreme condition. The muscle starts remodelling its inner structure by multiplication of fibers and their growth. The High Intensity Focused Electro-Magnetic uses a specific range of frequencies that does not allow muscle relaxation between two consecutive stimulations. The muscle is forced to remain in contracted state for multiple seconds. When repeatedly exposed to these high load conditions the muscle tissue is stressed and is forced to adapt. Recent studies reported that on average 15%- 16% increase in abdominal muscle thickness was observed in treated patients one to two months after High Intensity Focused Electro-Magnetic treatments.




H Type

Folding Type

Fat Type

After Childbirth



Build abdominal muscles
Perfect Firm abs
and V line abs



The EMS machine works on the abdomen area contracting and tightening the muscles to both burn fat and sculpt the area. This treatment can help you improve the muscle definition that you are already seeing from your fitness routine. It is also ideal for addressing that stubborn inch of fat on your stomach that just won't go away.



FOR THE ABDOMEN



EMS machine works on the buttocks as well, using the same supramaximal muscle contractions as with the abdomen area. Many patients opt to try this treatment in place of an invasive Brazilian butt lift. High Intensity Focused Electro-Magnetic technology delivers the world's first noninvasive butt lift procedure. This treatment can give you the curvature, tone, and lift that you desire. It is an alternative to the Brazilian Booty Lift procedure and provides a defined shape instead of adding excess fat to the area.



FOR THE BUTTOCK





Dedicated multi-directional handle fixed arm

Multi-directional fixed support arm can be adjusted freely at any angle of 360°

The multi-directional fixed arm greatly facilitates the treatment of various parts of the body, and perfectly solves the problem that the handles of the shoulders, handles, shoulder blades and other parts are not easy to fix. Users can easily enjoy the treatment of muscle gain and fat reduction in various parts of the body.



Intelligent Operating System

GENDER SELECTION


Male | Female

TREATMENT AREA SELECTION

Buttock | Abdomen | Thigh | Arm | Calf

WORKING MODE SELECTION

Hiit | Strength | Muscle Building | Shape Up | Fat Burning



The new intelligent operating system can automatically match different treatment parameters according to different treatment parts; at the same time, it is equipped with five different working modes to directly match the corresponding treatment parameters according to the treatment needs of customers. Extremely easy to operate

Comparison

MUSCULPTING NEO

MUSCULPTING PRO

Technology

High Intensity Focused Electro-Magnetic + RF

High Intensity Focused Electro-Magnetic

What can it do

Build & tone muscle , Burn Fat, Tighten Skin

Build & Tone Muscle

What does it treat

Muscle + Fat + Skin

Muscle

How long is the treatment

30minutes

30 minutes

How many sessions

4 minimum

4 minimum

Contactless or not

Needs to be on the skin tomake contact

Can be on the skin or overclothes

Discomfort

None

None

Recovery

Absolutely none

Absolutelly none

How often should one get it

Between 1-3 months, variable

Between 1-3 months, variable

BMI requirement

Up to BMI of 35

Up to BMI of 30

CONTACT US

E-mail: info@bjnubway.com

Phone: 0086-10-60418533

Web: www.nubway.com

Add: West zone building No.5 Mauhwa Industry Park No.1 Caida 3rd

Street Shunyi District Beijing 101300 China

SPECIFICATIONS

Frequency:	1-150hz
RF Power:	800W
Power:	5000w
Intensity:	7Tesla
Screen:	15inch tough screen
Pulse width:	310 μ s \pm 10%
Cooling:	Air Cooling
Voltage:	AC220V \pm 10% 10A 50HZ , 110v \pm 10% 10A 60HZ